

# Dual MassSTAR

LC-MS Sample Preparation with accelerated Performance



## Maximum Functionalities

Protein precipitation and solid-phase extraction on one system



## Standardized workflows


Verified protocols by Hamilton and customers



## Fast & Flexible

Processing of up to 80 samples below 1.5 h

# Specifications



- 1 Tip Carrier
- 2 Plate Carrier
- 3 Reagent Carrier
- 4 Sample Carrier
- 5 Hamilton Cooling Module
- 6 Hamilton Heater Shaker
- 7 CVS-Vacuum Station (optional)
- 8 Barcode Reader
- 9 Hamilton SBS Centrifuge

## Configurations and Ordering Information

Description	PN
Dual MassSTAR	806614
Options	
1. Immunosuppressants	
Roche Cobas Carrier	185037
2. Solid Phase Extraction	
CVS-Basic vacuum system 230V	199020
CVS waste bottle 4L	281540
DWP kit for CVS	190034
Insert kit MTP for CVS	190035

## Physical Dimensions

Description	
Width	1391 mm
Depth (w/o Autoload)	788 mm
Depth (with Autoload)	1008 mm
Height	901 mm

## Methods/Analytes

Description	Input Sample
Antiepileptics	Serum / EDTA-Plasma
Antibiotics	Serum / EDTA-Plasma
Antidepressants	Serum / EDTA-Plasma
Neuroleptics	Serum / EDTA-Plasma
Synthetic Cannabinoids	Urine
Immunosuppressants	Whole Blood

## Operating Conditions

Description	
Temperature	+15 - +35°C
Humidity	15 - 85 % w/o condensation
Altitude above sea level	0 - 2000 m

© 2023 Hamilton Company. All rights reserved. All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries. Lit. No. F-2110-05 – 09/2023

**HAMILTON**

To find a representative in your area, please visit:

[www.hamiltoncompany.com/contact](http://www.hamiltoncompany.com/contact)  
Email: [infoservice@hamiltonrobotics.com](mailto:infoservice@hamiltonrobotics.com)

**United States**  
+1-775-858-3000  
**United Kingdom, Ireland**  
+44 121 272 92 80  
**Brazil**  
+55 11 95914 5000  
**China**  
+86 21 6164 6567

**Denmark, Norway, Sweden, Finland**  
+46 8410 27 373  
**France**  
+33 184 008 420  
**Germany, Austria, Switzerland**  
+49 89 248 804 804

**Benelux**  
+31 40 209 178 0  
**Italy**  
+39 039 930 06 06  
**Japan**  
+81 3 6435 6850  
**Spain, Portugal**  
+34 930 186 262

