

Sample **Loading** to PCR Result **Analysis**

ExiStation™ FA 96/384



Next Generation Fully-Automated MDx System

“De-capping
is not your job!”

ExiStation™ FA 96/384

Hands-free!

All you do is insert the sample and it will do the rest



Auto De-capping Sample Loading

- Automated de-capping system for primary tubes which eliminates the risk of cross-contamination
- Automated sample preparation

Nucleic Acid Extraction & PCR setup

- Consistent and accurate DNA purification
- Automated PCR setup

Amplification & Detection

- Real-Time data analysis

BIONEER
Innovation • Value • Discovery

BIONEER Corporation
Bioneer Global Center, 71, Techno 2-ro,
Yuseong-gu, Daejeon, 34013, Republic of Korea
Tel: +82-42-939-6333
Fax: +82-42-939-6444
E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216 Oakland,
CA 94607, USA
Toll Free: +1-877-264-4300
Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

European Branch
Ludwig-Erhard-STR. 30-34, 65760
Eschborn, Germany
Tel: +49-6196-9699102
E-mail: euinfo@bioneer.com