



# This FACS makes the Difference

## What do you expect from your next clinical flow cytometer?

- ✓ **Quiet - small - automated** startup and shutdown?
- ✓ Only **one software for all assays** types - flexible and CE-IVD?
- ✓ **One click daily QC** checking >70 parameter for **reliable results**?
- ✓ **Bi-directional data transfer** - including PDF and for **all assay types**?
- ✓ **Unlimited amount of events** in a FCS file - acquired and analyzed?
- ✓ **Instrument setup** automatically **traces compensation**?
- ✓ **Automated pipetting steps** including tracking of all steps?
- ✓ Outlook to **fully automated sample to answer solution**?
- ✓ Microsoft® Windows® 10 LTSC IoT **improved data security**?




[bdbiosciences.com/en-eu](https://bdbiosciences.com/en-eu)

The BD FACSLytic™ flow cytometer and BD FACSDuet™ Sample Preparation System are for In Vitro Diagnostic Use. CE marked in compliance with the European In Vitro Medical Device Directive 98/79/EC  
BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2020 BD. All rights reserved.



# BD FACSLyric™ (CE-IVD)

Accuracy - Reproducibility - Simplicity





BD FACSLyric™ Flow Cytometer



BD OneFlow™ Solution



BD FACSTM Universal Loader



Flexible bidirectional data transfer



BD FACSLyric™ BD FACSDuet™ Sample Preparation System



BD Application Support

## Your advantages at a glance



### Consistent and reliable results

High sensitivity and resolution to resolve dim markers

Stable instrument performance with CV on <0.4%



### Small - compact - quiet - flexible

Find a place in every lab

Quiet operations with vacuum fluidics

Various data export formats to be used in alternative software



### Rare events at highest resolution

Outstanding resolution at all flow rates

No limit on events acquired

Minimal carry over to recover rare events



### Simplified operations

Pre-program startup and shutdown

One BD FACSuite™ software from daily QC, acquisition and analysis

CE-IVD matched reagent kits and software



### Standardization and reproducibility build in

Reproducible results over time and between instruments

Assay specific QC protocols & tracking

CE-IVD & FDA 21 CFR Part 11 support



### Ready for automation and bidirectional LIS connectivity

Universal tube and plate loader

Bi-directional LIS connectivity for all types of assays (CE-IVD and LDT)

BD FACSDuet Sample preparation

## We are your partner with experience



Your local BD Biosciences representative and our team of experts are committed to providing continuous support for your seamless laboratory routine.



### BD Customer Support Team

Telephone · E-Mail · Chat · Remote Support on demand

Detailed knowledge from sample preparation to data analysis



### On-site application support

Experienced BD application specialists provide individual on-site support in the areas of instruments, software and reagents.

## What customers are saying about our excellent support



"The installation of the new BD FACSLyric / FACSDuet system was supported by the local BD application specialist in an absolutely professional manner."

"All assays to be transferred and new ones to be prepared were established in close coordination and available for routine use immediately after the system was put into operation."



"Our Lyric has started and it all looks very good. Thank you very much for your perfect support during the transition."

"Reliable, always friendly and helpful guide through the variety of colors and extensive capabilities of the BD FACSLyric."

You like to convince yourself personally about the BD Support with our users? Please contact your BD representative and we will gladly connect you.